


EYECRYLTM PLUS

NATURAL YELLOW ASPHERIC HYDROPHOBIC ACRYLIC
FOLDABLE IOL WITH DELIVERY SYSTEM




FEATURES

- OPTIMALLY HYDRATED TO ENSURE "NO GLISTENINGS"
- NATURAL CHROMOPHORE
- 360 DEGREE ADVANCED SQUARE EDGE
- NEGATIVE S.A. TO COMPENSATE POSITIVE S.A. OF CORNEA
- EXCELLENT FOLDING & UNFOLDING PROPERTIES
- CONSISTENT BEHAVIOR IN ALL ENVIRONMENTAL CONDITIONS



SPECIFICATIONS

Model: ASHFY600



MATERIAL	Hydrophobic Acrylic containing Natural Chromophore		
OPTIC TYPE	Single Piece, 360° Square Edge with Aspheric Optic		
OPTIC SIZE	6.00 mm		
OVERALL SIZE	13.00 mm		
ANGULATION	0°		
ACD	5.28		
REFRACTIVE INDEX	1.48		
ESTIMATED A-CONSTANT	118.5		
RECOMMENDED OPTICAL A-CONSTANTS	SRK – T 118.6	SRK – II 118.9	Holl 1 Const SF : 1.63
	HOFFER Q ACD : 5.43		HAIGIS α0:1.20, α1:0.40, α2:1.10
	Barrett 1.67	Hill-RBF 118.60	Olsen ACD 4.47
DIOPTER RANGE	+5.0 to +30.0 D (with 0.5 D steps)		
IMPLANTATION SITE	Capsular Bag		
STERILIZATION	Irradiation		
SHELF LIFE	4 years from date of manufacture		

biotech
VISION CARE